# **E5CN**

**NAUDOJIMO INSTRUKCIJA  
INSTRUKCJA OBSŁUGI  
ИНСТРУКЦИЯ ПО ЭКСПЛУАТАЦИИ  
INSTRUCTION MANUAL  
BEDIENUNGSANLEITUNG**

**MRC LTD -HAGAVISH 3, HOLON ISRAEL**  
Tel: 972-3-5595252, Fax: 972-3-5594529  
E-mail: [mrc@mrclab.com](mailto:mrc@mrclab.com) • Web. Site: [www.mrclab.com](http://www.mrclab.com)

## TURINYS

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
## 1. Parts of the device and functions of the keys

PV No.1 display the present temperature in the chamber of furnace or a symbol.

SV No.2 display the set temperature, a symbol or value of parameter.

**OUT1** indicator shines when the heating elements are heating.

**STOP** indicator shines when the regulator is stopped.

Level key. This key switches on or off a level of adjustment.

Mode key. This key switches parameters.

Down key. Each pressing of this key reduces value of parameter or changes a symbol.

Up key. Each pressing of this key increases value of parameter or changes a symbol.

## 2. The control of device

The regulator has two levels of control: operating and adjustment.

Operating level is adjusted automatically at inclusion of device in a network. SV No.2 display info-line displays the set temperature, PV No.1 display info-line – present temperature in the chamber of furnace.

The level of regulation is used for setting of temperature control parameters.

When the voltage supply is off entered values of parameters remain.

## 3. Operating level

Parameters are switching by short pressing of parameters key. Values of parameters are adjusted by keys “down” ⏴ and “up” ⏵.

**SP-M**: displays temperature which should be in the chamber of furnace during set time. For example: if the temperature in chamber should be 200 °C **SP-M** displays: If speed

**SP-M**  
200

of heating is switched off at a level of regulation the device does not display this parameter.

**PRSE**: start of the program (status). If with “down ” or “up ” keys we will set: program will be stopped. This is a “reset” status.

**PRSE**  
**RSEE**

If we will set: **PRSE**  
**SESE** regulator will perform the set program and will switch off heating of furnace automatically. When the program is finished, in SV No.1 display the note “**END**” is blinking. The device is required to be set on the “reset” status.

**SKLR**: displays the time (in minutes) till finish of heating. For example: if up to the end of program there are 5 minutes the device displays **SKLR**  
**S** If temperature in the chamber of furnace


is not in the set range of temperature the certain time, this time is not included to the time of heating.

**RS**: heating switching on/off. If we will set: **R-S**  
**RUN** ( **R-S**  
**RSEL** ), furnace will heat up with the set speed to the fixed working temperature. Temperature will be constant until regulator will be switched off manually, set to **R-S**  
**STOP** status. Parameter **RL-2** value

should be equal to furnace maximum temperature

#### 4. Level of adjustment *L.AdJ*

Level of regulation is switching on/off by short pressing of the level key. Parameters are switching by short pressing of the parameters key. Values of parameters are setting by “down ” and “up ” keys.

**RL**: automatic selection of parameters (parameters of PID combination: **P**, **i**, **d**) of temperature regulation (further in text: AT). If the temperature in the chamber of furnace changes and is unstable it is required to switch on AT. AT can be switched on only when regulator executes the program.

When AT is switched: **RL**  
**RL-2** or **RL**  
**RL-1** SV No.2 display is blinking. The regulator switches off AT automatically when selects parameters of PID combination. It is recommended to select **P**, **i**, **d** parameters at operating temperature.

**ENS**: Movement of value of temperature which is shown on PV No 1display. It should be:

**ENS**  
**0.0**

$P$  : parameter of PID combination

$\bar{c}$  : parameter of PID combination

$d$  : parameter of PID combination

If parameters of PID combination are known, they can be entered manually or selected automatically (see AT). Parameters of PID combination depend on model of the furnace, operating temperature and loading.

$S\bar{a}Rk$  : time of keeping, in minutes (see diagram). For example: if the keeping time is 20 minutes, the device displays:  $\boxed{\begin{matrix} S\bar{a}Rk \\ 20 \end{matrix}}$ .

$Wt-b$  : range of temperature, °C. The range of temperature is related to time of heating (see diagram), but not related to accuracy of adjustment of temperature. It is recommended to set a range of temperature on 10 °C  $\boxed{\begin{matrix} Wt-b \\ 10 \end{matrix}}$ . In this case the time of keeping will be considered in

an interval of 10 degrees below and 10 degrees higher than the set temperature. At other values of temperatures the time is not included to the time of heating.

$SPRt$  : speed of heating, °C/minute. For example: if the speed of heating is 5 °C/min, the device  $\boxed{\begin{matrix} SPRt \\ 5 \end{matrix}}$

displays . If the speed of heating is switched off (the device displays  $\boxed{\begin{matrix} SPRt \\ off \end{matrix}}$ ), furnace is heating up with the fastest speed.

$\bar{a}L-H$  : the highest used power in percents. It should be 105.0 ( $\boxed{\begin{matrix} \bar{a}L-H \\ 105.0 \end{matrix}}$ ).

$\bar{a}L-L$  : the lowest used power in percents. It should be - 5.0 ( $\boxed{\begin{matrix} \bar{a}L-L \\ -5.0 \end{matrix}}$ ). If the value of  $\bar{a}L-L$  parameter will be positive the furnace will be heated up even at the switched off regulator.  $\bar{a}RL$  parameter value should be 0.0 ( $\boxed{\begin{matrix} \bar{a}RL \\ 0.0 \end{matrix}}$ ).

#### Notes:

1. If temperature in the chamber of furnace is 5 °C below set temperature (or the temperature which should be in a set time) the colour of PV No.1 display is yellow. If temperature in the chamber of furnace is 5 °C more than set temperature the colour of PV No.1 display is red. If the difference between present temperature and the set temperature is not more than 5 °C the colour of PV No.1 display is green.
2. If during 15 seconds any button will not be pressed, the device will automatically return to the initial status, it means: PV No.1 display displays present temperature in the furnace chamber, SV No.2 display – the set temperature.
3. Depending on furnace type some parameters could be not.

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
We will heat up a furnace up to 300 °C with a speed of 5 °C per minute, the reached temperature of 300±10°C will be constant within 80 minutes. After that heating of the furnace is switched off automatically.

**If during 15 seconds any key will not be pressed, the device will automatically return to the operating status.**

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
 <p>The diagram shows the Omron E5CN control panel. It features a central digital display showing '20' for PV (Present Value) and '0' for SV (Set Value). The display is flanked by '°C' on the left and 'PV' on the right for the top line, and 'SV' on the right for the bottom line. To the left of the display are four buttons: SUB1, SUB2, SUB3, and HA. Below these are two rows of buttons: the first row has OUT1, STOP, and a power button (a circle with a vertical line); the second row has OUT2, CMW, and MANU. At the bottom of the panel are four navigation buttons: a circle, a left arrow, a downward arrow, and an upward arrow. The model number 'E5CN' is printed at the bottom right of the panel.</p>	<p>1. PV infoeilutė rodo esamą temperatūrą kameroje 20°C SV infoeilutė užduotą temperatūrą 0°C 1. Wiersz informacyjny PV pokazuje obecną temperaturę w komorze 20°C SV infoeilutė užduotą temperatūrą 0°C 1. PV дисплей показывает настоящую температуру в камере 20°C SV дисплей показывает заданную температуру 0°C 1. PV info-line displays 20°C as a present temperature in the chamber. SV info-line displays 0°C as a set temperature. 1. Die Infozeile PV zeigt 20°C als die Ist-Temperatur im Heizraum an. Die Infozeile SV zeigt 0°C als die eingestellte Temperatur.</p>
--	--

**Reguliavimo lygmuo**  
**Poziom regulowania**  
**Режим настройки**  
**Level of regulation**  
**Regelungsebene**

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2. Short pressing of the level key  is switching on the level of adjustment.

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
3. Several times shortly press the parameters key  until the inscription **SōAK** will appear.

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
4. With the keys set the time of keeping.

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
5. Pressing of the key is switching on the settings of temperature range.

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
6. With the keys set the temperature range  $\pm 10,0^{\circ}\text{C}$ .

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
7. Pressing of the key is switching on the window of speed settings of heating up.

(If status of **SPRŁ** is **ōFF**, then furnace is heating up with the fastest speed)

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

8. With the keys we set 5 °C per minute as a speed of heating up.  
 PV info-line displays **SPRŁ**  
 SV info-line displays **5**

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
9. Pressing of the key is switching on the operating level.  
 PV info-line displays **20**  
 SV info-line displays **0**

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
10. With the   keys we set 300 °C as temperature of operating.

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