

**TRO-1200D, Digital Tube Roller, 5 Rollers**

Digital Roller Mixer TRO-1200D is your ideal choice for blending analysis samples such as blood samples. Gently rolling and tilting action prevents the sediment of the specimen to provide the homogenous specimen quality. Plastic rollers are suitable for use with a variety of tubes. Roller mixer TRO-1200D has a widely range of application.



TRO-1200D

**Features:**

- Brushless DC motor is durable and maintenance free.
- Precision timer can be set to 100 hours with 1 minute increment or timer can be set as a continuous mode.
- Digital speed setting can be set with 1rpm increment to 80 rpm.
- Gently rolling and tilting action.
- Easily load and unload without tuning of the unit.
- Plastic rollers are easy to clean.
- Quality certification: ISO 9001, ISO 13485.

Model	TRO-1200D
Speed	2~80rpm
Time Mode	Continuous mode / timed mode (100 hour timer)
Number of Rollers	5
Roller Dimensions	ø30 x L320mm
Overall Dimensions	W455 x D180 x H95 (mm)
Weight (N.W./G.W.)	3.7 / 4.7 Kgs
Power Supply	AC 100~240V 50/60Hz