



MX-S

### MX-S, Adjustable speed Tube Vortex

- Touch operation or continuous mode
- Variable speed control from 0 to 3000rpm
- Used for various mixing applications with optional adapters
- Specially designed vacuum suction feet for body stability
- Robust aluminum-cast construction.

Model	MX-S
<b>Voltage</b>	110-120V/220-230V,50/60Hz
<b>Power</b>	60W
<b>Mixing motion</b>	Orbital
<b>Orbital diameter</b>	4mm
<b>Motor type</b>	Shaded pole motor
<b>Motor rating input / output</b>	58W / 10W
<b>Speed range</b>	100-3000rpm
<b>Speed display</b>	Scale
<b>Run type</b>	Touch operation/Continuous
<b>Dimension (W×H×D)</b>	127×130×160mm
<b>Weight</b>	3.5kg
<b>Permissible ambient temperature</b>	5-40°C
<b>Permissible relative humidity</b>	80%RH
<b>Protection class</b>	IP21