

# AUTOCLAVES

## Table Top



**STE-F-8**

### STE-F-8, 8 Liter Autoclave

#### Default Program:

Program	Temperature	Pressure	Holding time (min)
Quick 134°C	134°C	210 Kpa	3.5
Wrapped 134°C	134°C	210 Kpa	4
Wrapped 121°C	121°C	110 Kpa	20
Dry	-	-	-
B&D Test	134°C	210 Kpa	3.5
Helix Test	134°C	210 Kpa	3.5
Vacuum test	-	-	-

Model	Chamber Size	Type	Capacity (L)	Overall Size(WxHxD)	Packing Size (WxHxD)	Weight (kg)
STE-F-8	Ø170x320mm 2pcs, Max 3pcs (Tray's optional)	Class B	8L	420x370x525mm	503x438x640mm	37